

NUF2 antibody

Catalog Number: 113365

Product name

NUF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NUF2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, human liver tissue, human skeletal muscle tissue, human testis tissue, PC-3 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of NUF2 recombinant protein (Accession Number: XM_024450112). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113365(NUF2 antibody) at dilution of 1:500



IP Result of anti-NUF2 (IP:Catalog No:113365, 3ug; Detection:Catalog No:113365 1:500) with HEK-293 cells lysate 1760ug.