

# RNFT2 antibody



Catalog Number: 114765

## Product name

RNFT2 antibody

## Specificity

Human; other species not tested.

## Antibody description

RNFT2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, SMMC-7721 cells. Observed molecular weight by Western-blot: 49-59 kDa

## Preparation

This antibody was obtained by immunization of RNFT2 recombinant protein (Accession Number: BC011878). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114765(RNFT2 Antibody) at dilution of 1:300