

# Anti-RPS6 Rabbit antibody

Catalog Number: 166765

## Product name

Anti-RPS6 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

RPS6 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human RPS6

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/1000

IHC: 1/200


ICC/IF: 1/100

FC: 1/100

## Validations



Immunocytochemistry of RPS6(green) in HeLa cells using RPS6 Rabbit pAb at dilution 1/50, and



Western blot detection of RPS6 in K562,C6,3T3 cell lysates using RPS6 Rabbit pAb(1:1000 diluted).Predicted band size:29kDa.Observed band size:29kDa.