

ADAM2 antibody



Catalog Number: 107765

Product name

ADAM2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ADAM2 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human testis tissue. Observed molecular weight by Western-blot: 50-55 kDa, 65-75 kDa

Preparation

This antibody was obtained by immunization of ADAM2 recombinant protein (Accession Number: NM_001278113). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107765(ADAM2 antibody) at dilution of 1:1000