

Anti-NFIB Rabbit antibody



Catalog Number: 168765

Product name

Anti-NFIB Rabbit antibody

Specificity

Human, Rat

Antibody description

NFIB Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human NFIB

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of NFIB in C6,HeLa cell lysates using NFIB Rabbit mAb(1:1000 diluted).Predicted band size:47kDa.Observed band size:47kDa.