

DCUN1D5 antibody



Catalog Number: 109765

Product name

DCUN1D5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCUN1D5 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, mouse brain tissue, rat spleen tissue, U-937 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of DCUN1D5 recombinant protein (Accession Number: NM_032299). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109765(DCUN1D5 antibody) at dilution of



1:500

IP Result of anti-DCUN1D5 (IP:Catalog No:109765, 3ug; Detection:Catalog No:109765 1:800) with HeLa cells lysate 1840ug.