

DMXL1 antibody

Catalog Number: 110075

Product name

DMXL1 antibody

Specificity

Human; other species not tested.

Antibody description

DMXL1 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells.

Preparation

This antibody was obtained by immunization of DMXL1 recombinant protein (Accession Number: NM_001349239). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-DMXL1 (IP:Catalog No:110075, 4ug; Detection:Catalog No:110075 1:300) with HEK-293 cells lysate 1480ug.