

KLC3 antibody



Catalog Number: 112075

Product name

KLC3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KLC3 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of KLC3 recombinant protein (Accession Number: BC025318). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112075(KLC3 antibody) at dilution of 1:500