

Anti-MAP2K6 (8A11) Mouse antibody



Catalog Number: 168075

Product name

WB, IP

Anti-MAP2K6 (8A11) Mouse antibody

Dilutions

Specificity

WB: 1/1000

Human, Rat

Antibody description

MAP2K6 (8A11) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human MKK6 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Western blot analysis of extracts from PC12, C6, Ramos, Jurkat and HeLa cell lysates using MKK3/6 mouse mAb (1:1000 diluted). Predicted band size: 40KDa. Observed band size: 40KDa.



Immunoprecipitation analysis of HeLa cell lysates using MKK3/6 mouse mAb.