

# Anti-FOXO1 Rabbit antibody



Catalog Number: 169075

## Product name

Anti-FOXO1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FOXO1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human FoxO1a

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/1000-1/2000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200


## Validations



Immunohistochemical analysis of paraffin-embedded human brain carcinoma, using FoxO1a Antibody.



Immunofluorescent analysis of HeLa cells, using



FoxO1a Antibody .

Western blot detection of FoxO1a in Jurkat, THP-1, SKOV-3, MCF7 and HeLa cell lysates using FoxO1a rabbit pAb (1:1000 diluted). Predicted band size: 70kDa. Observed band size: 78 to 82kDa.