

SLC26A11 Polyclonal Antibody



Catalog Number: 163475

Product name

SLC26A11 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to SLC26A11

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 500-600 of human SLC26A11 (NP_001159819.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - SLC26A11 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SLC26A11 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.