

Phospho-MKKNK1-T197/202 pAb



Catalog Number: 166475

Product name

Phospho-MKKNK1-T197/202 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-MKKNK1-T197/202

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding T197/202 of human MKKNK1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-MKKNK1-T197/202 pAb

Western blot analysis of extracts of various cell lines, using Phospho-MKKNK1-T197/202 antibody at 1:1000 dilution. HeLa cells were treated by 10% FBS for after serum-starvation overnight or treated by ATP (5mM) at 30u2103 for 1 hour .NIH/3T3 cells were treated by 10% FBS for 30 minutes after serum-starvation overnight or treated by ATP (5mM) at 30u2103 for 1 hour .PC-12 cells were treated by 20% FBS for 20 minutes after serum-starvation overnight or treated by ATP (5mM) at 30u2103 for 1 hourSecondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.Detection: ECL Basic Kit .Exposure time: 1s.