

Anti-LC3B Rabbit antibody



Catalog Number: 167475

Product name

Anti-LC3B Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

LC3B Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human LC3B

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/500-1/2000


ICC/IF: 1/50-1/200

IP: 1/50

Validations



Immunofluorescent analysis of HeLa cells treated with chloroquine, using LC3B Antibody.



Western blot analysis of LC3B expression in (1) Human brain lysate; (2) RAW 264.7 cell lysate; (3) C6 cell lysate.