

BNP antibody

Catalog Number: 108475

Product name

BNP antibody

Specificity

Human; other species not tested.

Antibody description

BNP Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human heart tissue.

Preparation

This antibody was obtained by immunization of BNP recombinant protein (Accession Number: NM_002521). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:117210(BNP Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:117210(BNP Antibody) at dilution of 1:50