

# FOSL2 Polyclonal Antibody



Catalog Number: 163875

## Product name

FOSL2 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to FOSL2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 187-326 of human FOSL2 (NP\_005244.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200


IF 1:50 - 1:200

## Validations



## Western blot - FOSL2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using FOSL2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



## Immunohistochemistry - FOSL2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human esophagus using FOSL2 antibody at dilution of 1:100 (40x lens).

# FOSL2 Polyclonal Antibody



Catalog Number: 163875

---



Immunofluorescence - FOSL2 Polyclonal Antibody

Immunofluorescence analysis of U-2 OS cells using FOSL2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.