

# NCR1 Polyclonal Antibody



Catalog Number: 164875

---

## Product name

NCR1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to NCR1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-258 of human NCR1 (NP\_001138929.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

## Validations



Immunofluorescence - NCR1 Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using NCR1 antibody . Blue: DAPI for nuclear staining.