

TBL1X antibody

Catalog Number: 115875

Product name

TBL1X antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TBL1X Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue. Observed molecular weight by Western-blot: 57kd

Preparation

This antibody was obtained by immunization of TBL1X recombinant protein (Accession Number: NM_001139468). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:115875(TBL1X antibody) at dilution of 1:400