

# Anti-hnRNP G-T Rabbit antibody



Catalog Number: 168875

## Product name

Anti-hnRNP G-T Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

hnRNP G-T Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human hnRNP G-T

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of hnRNP G-T in K562, Rat Brain, C6, 3T3, HeLa cell lysates using hnRNP G-T Rabbit mAb(1:1000 diluted). Predicted band size: 43kDa. Observed band size: 43kDa.