

Anti-hnRNP G-T Rabbit antibody



Catalog Number: 168875

Product name

Anti-hnRNP G-T Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

hnRNP G-T Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human hnRNP G-T

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of hnRNP G-T in K562, Rat Brain, C6, 3T3, Hela cell lysates using hnRNP G-T Rabbit mAb(1:1000 diluted). Predicted band size: 43kDa. Observed band size: 43kDa.