

EIF1AX antibody



Catalog Number: 110175

Product name

EIF1AX antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EIF1AX Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A375 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 16-23 kDa

Preparation

This antibody was obtained by immunization of EIF1AX recombinant protein (Accession Number: XM_016999709). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110175(EIF1AX Antibody) at dilution of 1:600



Immunofluorescent analysis of HepG2 cells using Catalog No:110175(EIF1AX Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)