

LCN9 antibody

Catalog Number: 112175

Product name

LCN9 antibody

Specificity

Human; other species not tested.

Antibody description

LCN9 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of LCN9 recombinant protein (Accession Number: NM_001001676). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112175(LCN9 Antibody) at dilution of 1:700