

Anti-EphB6/Eph Receptor B6 antibody



Catalog Number: 105175

Product name

Anti-EphB6/Eph Receptor B6 antibody

Specificity

Human EPHB6

Antibody description

Rabbit polyclonal to EphB6/Eph Receptor B6

Preparation

Produced in rabbits immunized with A synthetic peptide corresponding to the N-terminus of the human EPHB6., and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 0.5-1 µg/mL

Validations



Lanes	A
Items	
Sample (Recombinant protein)	EPHB6
Sample Volume (ng/lane)	100
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	0.5-1µg/ml
Secondary Antibody	HRP Conjugated Goat anti-rabbit IgG (H+L) at 1:1000 dilution

EPHB antibody, Rabbit PAb, Antigen Affinity Purified, Western blot