

# TMEM201 antibody



Catalog Number: 116175

## Product name

TMEM201 antibody

## Immunogen

[Human TMEM201 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TMEM201 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 43-45 kDa

## Preparation

This antibody was obtained by immunization of TMEM201 recombinant protein (Accession Number: XM\_011540910). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116175(TMEM201 Antibody) at dilution of 1:100