

# ZNF419 antibody



Catalog Number: 117175

---

## Product name

ZNF419 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZNF419 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, Jurkat cells. Observed molecular weight by Western-blot: 58 kDa

## Preparation

This antibody was obtained by immunization of ZNF419 recombinant protein (Accession Number: NM\_024691). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:117175(ZNF419 Antibody) at dilution of 1:600