

# Anti-H-FABP (5C4F5) Mouse antibody



Catalog Number: 168175

## Product name

Mouse IgG1

Anti-H-FABP (5C4F5) Mouse antibody

## Applications

## Specificity

sELISA Capture

Human

## Validations

## Antibody description

H-FABP (5C4F5) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant H-FABP protein expressed in E.coli


## Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Standard Curve for H-FABP: Capture Antibody Mouse mAb (5C4-F5) to H-FABP at 4u03bcg/ml and Detector Antibody Mouse mAb(3E8-B10u30017B10-G11u30014H1-B8)to H-FABP at 0.08 u03bcg/ml.