

CDT1 antibody



Catalog Number: 109175

Product name

CDT1 antibody

Specificity

Human; other species not tested.

Antibody description

CDT1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of CDT1 recombinant protein (Accession Number: NM_030928). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109175(CDT1 antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:109175(CDT1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)