

# PAPD5 antibody

Catalog Number: 113575

## Product name

PAPD5 antibody

## Specificity

Human; other species not tested.

## Antibody description

PAPD5 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 63-80 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001040285). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:2000

## Validations



IP Result of anti-PAPD5 (IP:Catalog No:113575, 3ug; Detection:Catalog No:113575 1:500) with Jurkat cells lysate 1600ug.