

APOC1 Polyclonal Antibody



Catalog Number: 164575

Product name

APOC1 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to APOC1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-83 of human APOC1 (NP_001636.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:100

Validations



Immunohistochemistry - APOC1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat liver using APOC1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - APOC1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using APOC1 antibody at dilution of 1:100 (40x lens).