

Anti-SIRP alpha Rabbit antibody



Catalog Number: 167575

Product name

Anti-SIRP alpha Rabbit antibody

Specificity

Human

Antibody description

SIRP alpha Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human SIRP alpha

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

FC: 1/50

Validations



Western blot detection of SIRP alpha in A549, MCF-7, 293T cell lysates using SIRP alpha Rabbit mAb(1:1000 diluted). Predicted band size: 55kDa. Observed band size: 85kDa.