

# Anti-FTO Rabbit antibody



Catalog Number: 168575

## Product name

Anti-FTO Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IHC-P

## Antibody description

FTO Rabbit monoclonal antibody

## Dilutions

1:1000

## Preparation


Antigen: A synthetic peptide of human FTO

1:20

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Validations



## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of FTO in K562, HeLa cell lysates using FTO Rabbit mAb(1:1000 diluted). Predicted band size: 58kDa. Observed band size: 58kDa.

## Clonality

Monoclonal