

Human APOA1BP / AIBP (His Tag) recombinant protein



Catalog Number: 501975

General Information

Gene Name Synonym

Apolipoprotein A-I-binding protein; NAD(P)HX epimerase; YjeF N-terminal domain-containing protein 1

Protein Construction

A DNA sequence encoding the human APOA1BP (CAC86580.1) (Met1-Gln288) was expressed with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 25

Molecular Mass

The recombinant human APOA1BP consists 275 amino acids and predicts a molecular mass of 30.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human APOA1BP / AIBP Protein (His Tag)