

DOCK10 antibody

Catalog Number: 109975

Product name

DOCK10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DOCK10 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 250 kDa

Preparation

This antibody was obtained by immunization of DOCK10 recombinant protein (Accession Number: BC015018). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-Dock10 (IP:Catalog No:109975, 4ug; Detection:Catalog No:109975 1:300) with HeLa cells lysate 3200ug.