

Human IL17B/ZCYTO7 (Fc Tag) recombinant protein



Catalog Number: 500275

General Information

Gene Name Synonym

Cytokine Zcyto7; Interleukin-20; Neuronal interleukin-17-related factor

Protein Construction

A DNA sequence encoding the human IL17B (NP_055258.1) (Met1-Phe180) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 21

Molecular Mass

The recombinant human IL17B consists 398 amino acids and predicts a molecular mass of 44.9 kDa.

Formulation

Lyophilized from sterile 100 mM Glycine, 10 mM NaCl, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IL17B / IL-17B Protein (Fc Tag)