

# SHKBP1 antibody



Catalog Number: 115275

## Product name

SHKBP1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

SHKBP1 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 76kd

## Preparation

This antibody was obtained by immunization of SHKBP1 recombinant protein (Accession Number: NM\_138392). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:115275(SHKBP1 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using SHKBP1 antibody Catalog No:115275 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).