

Anti-STAT5B (7H9) Mouse antibody



Catalog Number: 168275

Product name

Anti-STAT5B (7H9) Mouse antibody

Specificity

Human

Antibody description

STAT5B (7H9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human STAT5B expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, IP

Dilutions


WB: 1/500 - 1/2000

ELISA: 1/10000

Validations



Western blot detection of stat5b in Jurkat, PC-12, C6, 3T3 and HeLa cell lysates using stat5b mouse mAb (1:1000 diluted). Predicted band size: 90KDa. Observed band size: 90KDa.



Immunoprecipitation analysis of HeLa cell lysates using stat5b mouse mAb.