

IGBP1 antibody



Catalog Number: 111675

Product name

IGBP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

IGBP1 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HeLa cells, Jurkat cells, MCF7 cells. Observed molecular weight by Western-blot: 40-45kd

Preparation

This antibody was obtained by immunization of IGBP1 recombinant protein (Accession Number: NM_001551). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111675(IGBP1 antibody) at dilution of 1:400