

Recoverin antibody



Catalog Number: 114675

Product name

Recoverin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Recoverin Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue, Y79 cells. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of Recoverin recombinant protein (Accession Number: NM_002903). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:114675(RCVRN antibody) at dilution of 1:400



Immunofluorescent analysis of MCF-7 cells, using RCVRN antibody Catalog No:114675 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).