

# Stanniocalcin 2 antibody



Catalog Number: 115675

## Product name

Stanniocalcin 2 antibody

## Specificity

Human; other species not tested.

## Antibody description

Stanniocalcin 2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human breast cancer tissue, human kidney tissue, human renal cell carcinoma tissue, human stomach cancer tissue. Observed molecular weight by Western-blot: 33 kDa

## Preparation

This antibody was obtained by immunization of Stanniocalcin 2 recombinant protein (Accession Number: NM\_003714). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115675(STC2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115675(STC2 antibody) at dilution of 1:100 (under 10x lens)