

# Anti-IGHA1 antibody



Catalog Number: 176675

## Product name

Anti-IGHA1 antibody

## Specificity

Human

## Antibody description

Mouse monoclonal antibody to IGH A1

## Preparation

This antigen of this antibody was native protein.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

Elisa, WB

## Dilutions

Elisa:1:5000

WB:1:500-1:2,000

## Validations




Fig1: Western blot analysis of IgA on human plasma lysate. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody was used at a 1:500 dilution in 5% BSA at room temperature for 2 hours. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1:5,000 dilution was used for 1 hour at room temperature.