

C1orf26 antibody

Catalog Number: 108675

Product name

C1orf26 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

C1orf26 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 100-105 kDa

Preparation

This antibody was obtained by immunization of C1orf26 recombinant protein (Accession Number: NM_017673). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:2000

Validations



IP Result of anti-C1orf26 (IP:Catalog No:108675, 5ug; Detection:Catalog No:108675 1:500) with mouse testis tissue lysate 8000ug.