

Cyclophilin E antibody



Catalog Number: 109675

Product name

Cyclophilin E antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Cyclophilin E Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of Cyclophilin E recombinant protein (Accession Number: NM_006112). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109675(PPIE Antibody) at dilution of 1:1500