

# NUP153 antibody

Catalog Number: 113375

## Product name

NUP153 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NUP153 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, A431 cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 153 kDa

## Preparation

This antibody was obtained by immunization of NUP153 recombinant protein (Accession Number: NM\_005124). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113375(NUP153 antibody) at dilution of 1:500



IP Result of anti-NUP153 (IP:Catalog No:113375, 4ug; Detection:Catalog No:113375 1:500) with K-562 cells lysate 3600ug.