

SIGLEC7 Polyclonal Antibody



Catalog Number: 164375

Product name

SIGLEC7 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SIGLEC7

Preparation

Antigen: Recombinant protein of human SIGLEC7

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:20 - 1:200

Validations



Western blot - SIGLEC7 Polyclonal Antibody

Western blot analysis of extracts of A-549 cells, using SIGLEC7 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.