

ROBO3-Specific antibody

Catalog Number: 114775

Product name

ROBO3-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ROBO3-Specific Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue, human brain tissue. Positive WB detected in mouse brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_022370). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114775(ROBO3-Specific antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human gliomas slide using Catalog No:114775(ROBO3-Specific Antibody) at dilution of 1:50

