

# Anti-Acetyl-Histone H3 (Lys9) Rabbit antibody



Catalog Number: 166775

## Product name

Anti-Acetyl-Histone H3 (Lys9) Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Acetyl-Histone H3 (Lys9) Rabbit polyclonal antibody

## Preparation

Antigen: The antiserum was produced against synthesized peptide derived from human Histone H3 around the acetylated site of Lys9. AA range:3-52

## Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, ELISA


## Dilutions

WB: 1/500-2000

IHC-p: 1/50-300

IF: 1/50-300

## Validations



Western blot detection of Acetyl-Histone H3 (Lys9) in COS7 and C6 cell lysates using Acetyl-Histone H3 (Lys9) rabbit pAb (1:3000 diluted). Predicted band size:17kDa. Observed band size:17kDa.