

# ADARB1 antibody



Catalog Number: 107775

## Product name

ADARB1 antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

ADARB1 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, A431 cells. Positive IP detected in A431 cells. Observed molecular weight by Western-blot: 70 kDa, 80 kDa

## Preparation

This antibody was obtained by immunization of ADARB1 recombinant protein (Accession Number: NM\_001112). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107775(ADARB1 antibody) at dilution of



IP Result of anti-ADARB1 (IP:Catalog No:107775, 4ug; Detection:Catalog No:107775 1:1000) with A431 cells lysate 1200ug.