

PPARGC1B antibody



Catalog Number: 114085

Product name

PPARGC1B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PPARGC1B Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, Transfected HEK-293 cells. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 100-110 kDa

Preparation

This antibody was obtained by immunization of PPARGC1B recombinant protein (Accession Number: NM_133263). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114085(PPARGC1B Antibody) at dilution of



1:300

IP Result of anti-PPARGC1B (IP:Catalog No:114085, 4ug; Detection:Catalog No:114085 1:300) with mouse liver tissue lysate 4000ug.