

TMED2 antibody

Catalog Number: 116085

Product name

TMED2 antibody

IHC: 1:20-1:200


Specificity

Human, Mouse; other species not tested.

Antibody description

TMED2 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in mouse pancreas tissue, mouse liver tissue, mouse skin tissue. Observed molecular weight by Western-blot: 23 kDa

Validations




Preparation

This antibody was obtained by immunization of TMED2 recombinant protein (Accession Number: NM_006815). Purification method: Antigen affinity purified.

mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:116085(TMED2 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:116085(TMED2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:116085(TMED2 Antibody) at dilution of 1:50

TMED2 antibody



Catalog Number: 116085

(under 40x lens)