

Anti-EpCAM (1D5) Mouse antibody



Catalog Number: 168085

Product name

Anti-EpCAM (1D5) Mouse antibody

Specificity

Human

Antibody description

EpCAM (1D5) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant EpCAM protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB, IP


Dilutions

WB: 1/1000

Validations



Immunoprecipitation analysis of HCT116 and HeLa cell lysates using EpCAM.



Western blot analysis of extracts from CHO-K1, CHO-K1 transfected by Flag-EpCAM, MCF7(EpCAM positive), HCT116(EpCAM positive), and HeLa(EpCAM negative) cell lysates using EpCAM mouse mAb (1:1000 diluted). Predicted band size: 39KDa. Observed band size: 39/35KDa.