

# Anti-EpCAM (1D5) Mouse antibody



Catalog Number: 168085

## Product name

Anti-EpCAM (1D5) Mouse antibody

## Dilutions

WB: 1/1000

## Specificity

Human

## Antibody description

EpCAM (1D5) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant EpCAM protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

WB, IP

## Validations



Immunoprecipitation analysis of HCT116 and HeLa cell lysates using EpCAM.



Western blot analysis of extracts from CHO-K1, CHO-K1 transfected by Flag-EpCAM, MCF7 (EpCAM positive), HCT116 (EpCAM positive), and HeLa (EpCAM negative) cell lysates using EpCAM mouse mAb (1:1000 diluted). Predicted band size: 39KDa. Observed band size: 39/35KDa.